

'Sup Sip-n-Puff user manual

Mouse actions:

Action:

Input:

Uses:

Hard sip	Right click	Options for folders, files or web links
Soft sip	Scroll down	To move a webpage down
Soft puff	Scroll up	To move a webpage up
Hard puff	Left click	Open programs, files, links, or items.
2 short hard puffs	Double click	To open programs, files or folders
Holding a hard puff	Clicking & dragging	To move items, files, folders, or programs

To keep your sip-n-puff safe and clean, we recommend running the mouthpiece through a dishwasher on a normal cycle at least once a week.

To disconnect the mouthpiece, just pull the tubing off the mouthpiece, then pull the mouthpiece directly off the device. Installation is the reverse of removal.

